Gelatin

Intended Use

Gelatin is used in preparing microbiological culture media.

Summary and Explanation

Gelatin is a protein of uniform molecular constitution derived chiefly by the hydrolysis of collagen.¹ Collagens are a class of albuminoids found abundantly in bones, skin, tendon, cartilage and similar animal tissues.1

Koch¹ introduced gelatin into bacteriology when he invented the gelatin tube method in 1875 and the plate method in 1881. This innovation, a solid culture method, became the foundation for investigation of the propagation of bacteria.¹ However, gelatinbased media were soon replaced by media containing agar as the solidifying agent.

Gelatin is used in culture media for determining gelatinolysis (elaboration of gelatinases) by bacteria. Levine and Carpenter² and Levine and Shaw³ employed gelatin media in their studies of gelatin liquefaction. Garner and Tillett⁴ used culture media prepared with gelatin to study the fibrinolytic activity of hemolytic streptococci.

Gelatin is a high grade gelatin in granular form which may be used as a solidifying agent or may be incorporated into culture media for various uses. Gelatin is used in Nutrient Gelatin, Motility GI Medium, Stock Culture Agar and Dextrose Starch Agar. A 0.4% gelatin medium is used in the presumptive differentiation of Nocardia brasiliensis from N. asteroides (see Nocardia Differentiation Media). Media containing gelatin are specified in standard methods^{5,6} for multiple applications.

Principles of the Procedure

The melting point of a 12% concentration of gelatin is between 28 and 30°C, which allows it to be used as a solidifying agent. Certain microorganisms elaborate gelatinolytic enzymes (gelatinases) which hydrolyze gelatin, causing liquefaction of a solidified medium or preventing the gelation of a medium containing gelatin. Gelatin is also used as a source of nitrogen and amino acids.

Procedure

See appropriate references for specific procedures using gelatin.

Expected Results

Refer to appropriate references and procedures for results.

References

- Gershenfeld and Tice. 1941. J. Bacteriol. 41:645. Levine and Carpenter. 1923. J. Bacteriol. 8:297.
- 3
- Levine and Shaw. 1924. J. Bacteriol. 9:225. Garner and Tillett. 1934. J. Exp. Med. 60:255.
- 5. U.S. Food and Drug Administration. 2001. Bacteriological analytical manual, online. AOAC International, Gaithersburg, Md.
- Eaton, Rice and Baird (ed.). 2005. Standard methods for the examination of water and wastewater, 21st ed., online. American Public Health Association, Washington, D.C. 6.

Availability

Difco[™] Gelatin

Cat. No.	214340	Dehydrated – 500 g
	214320	Dehydrated – 10 kg

User Quality Control

Identity Specifications

Difco™ Gelatin	
Dehydrated Appearance:	Light beige, free-flowing, homogeneous.
Solution:	12% solution, soluble in purified water upon slight heat- ing in a 50-55°C water bath. Solution is light amber, clear to slightly opalescent, may have a slight precipitate.
Prepared Gel:	Very light amber, clear to slightly opalescent, may have a slight precipitate.
Reaction of 2% Solution at 25°C:	рН 6.8 ± 0.2

Cultural Response Difco[™] Gelatin

Prepare a 12% Gelatin solution in 0.8% Nutrient Broth. Dispense into tubes and autoclave. Inoculate and incubate at 35 ± 2°C under appropriate atmospheric conditions for 18-48 hours or for up to 2 weeks for the gelatinase test. To read gelatinase, refrigerate until well-chilled and compare to uninoculated tubes. Tubes positive for gelatinase will remain liquid.

		0110 1110 101	
10 ² -10 ³	Good	+	
7 10 ² -10 ³	Good	+	
2 10 ² -10 ³	Good	-	
	$\begin{array}{c} 10^2 - 10^3 \\ 7 & 10^2 - 10^3 \\ 2 & 10^2 - 10^3 \end{array}$	10 ² -10 ³ Good 7 10 ² -10 ³ Good 2 10 ² -10 ³ Good	Cross Record (100 - 100) 10 ² -10 ³ Good + 7 10 ² -10 ³ Good + 2 10 ² -10 ³ Good -



